

# Low Vision Glasses Series

## Prismatic Reading Glasses



**Two Lens System:** in the two lens system, one lens is high powered “plus” lens for magnification, and the other a “prism” to alleviate the stress of converging the two eyes. The two lenses are bonded together. The eyeglass prescription can be incorporated into the system. The system was designed for tasks that require a working distance of between four and ten inches.

**One Lens Half Eye System:** the one piece prismatic system is designed as a lighter weight, more cosmetic alternative to the two lens system. It is intended for those with a minimal eyeglass prescription.

Some of the tasks and activities our patients have used their Prismatic Reading Glasses for include:

- Reading & writing
- Viewing photos
- Performing hand crafts
- Performing kitchen activities
- Polishing nails

# Low Vision Glasses Series

## Prismatic Reading Glasses



**“Very comfortable glasses! I use them for writing checks and doing my cross word puzzles.” Donald R.**

**“I can actually peel potatoes and see what I’m doing! I love it. Next I’ll try to file my nails.” Clara T.**



**“These glasses let me read the sports pages again. I can see to sign my name. Life is not as frustrating anymore.” Leo C.**